

9252120

## Wireless Gen2 2 Button Bed Station

### Description


The Wireless Gen2 2 Button Two Bed Station is a multifunction panel that can transmit alerts wirelessly to a miPOSITION enabled miALERT Nurse Call System using the proprietary miNET wireless transmission standard. Alerts can be initiated using the face button, or any device connected via quarter inch jack, such as a bed cord or bed pad. The device has a physical cancel button.

The unit mounts to a 1-gang box and the wireless design allows for easy installation with no cable runs. Low battery reporting is available through the miALERT Event Management Software.

These panels can also be used with Gen2 wireless enabled miALERT voice to room devices.

### Features

- Physical cancel from the unit.
- 1/4" jack with cord-out alarming.
- Wireless transmission requires no physical cabling.
- Battery powered, no direct power required.
- Two year battery life under nominal conditions.



### Specifications

**Unit:** 105mm x 53mm x 49mm

**Unit with mounting plate:** 115mm x 70mm x 49mm

**Unit exposure from wall:** 115mm x 70mm x 12mm

**Power:** Single 3V CR2 battery cell. Expected battery life under normal conditions is approximately 2 years. Units should be tested at least once a month to ensure battery condition is nominal. For low battery reporting a miALERT Event Management Software server must be present.

Product is designed to conform to the UL 2560 standard for wireless nurse call devices.

### Package Contents

- miBUS 2-Button Bed Station
- Faceplate
- Two (2) 6/32x1 Gang Box Screws
- One (1) CR2 3V Lithium Battery