

9272020

Wireless Gen2 2 Button Pull Cord

Description


The Wireless Gen2 Pull Cord is a reliable pull cord solution that can transmit alerts wirelessly to a miPOSITION enabled miALERT Nurse Call System using the proprietary miNET wireless transmission standard. Alerts are initiated using the pull cord and can be cancelled with a physical button the unit face.

The unit mounts to a 1-gang box and the wireless design allows for easy installation with no cable runs. Alerts can be transmitted up to 100ft under nominal conditions. Low battery reporting is available through the miALERT Event Management Software.

These pull cords can also be used with Gen2 wireless enabled miALERT voice to room devices.

Features

- Physical cancel from the unit.
- Wireless transmission requires no physical cabling.
- Battery powered, no direct power required.
- Two year battery life under nominal conditions.



Specifications

Unit: 105mm x 53mm x 36mm

Unit with mounting plate: 115mm x 70mm x 36mm

Unit exposure from wall: 115mm x 70mm x 15mm

Power: Single 3V CR2 battery cell. Expected battery life under normal conditions is approximately 2 years. Units should be tested at least once a month to ensure battery condition is nominal. For low battery reporting a miALERT Event Management Software server must be present.

Product is designed to conform to the UL 2560 standard for wireless nurse call devices.

For replacement batteries, contact NSI Distribution (www.nsidistribution.com).

Package Contents

- miBUS 4-Button Bed Station
- Faceplate
- Two (2) 6/32x1 Gang Box Screws
- One (1) CR2 3V Lithium Battery